


WINLAB PTY LTD

2 Pinnacle St, Brendale, Queensland, Australia, 4500
PO Box 5007, Brendale, Queensland, Australia, 4500
Ph: +61 7 3205 1209
Email: info@winlab.com.au www.winlab.com.au




Histology – Cytology – Adaptor Staining Hangers

OEM original and replacement hangers and Carriers integrate many different Slide Racks to a large range of staining machines:

- Dako®
- Hacker®
- InstrumeC®
- Leica®
- Medite® - Cot 20 & TST44
- Philips Digital Pathology®
- Raymond Lamb®
- Sakura®
- Stainmate®
- Tharmac®
- Thermo Fisher Scientific®
- Triangle Biomedical®
- Ventana®
- Vision BioSystems® and more

For more information, and your freight inclusive quote today, contact info@winlab.com.au



AWS101-3-1: Adaptor Hanger, carries all **Winlab® SR30 Series Racks**. Compatible with **Leica®, Thermo®, Tharmac®** and **Raymond Lamb®** staining machines.

AWS101-3-1D: Carries **InstrumeC®** and **Dako® 30 Place Racks**. Compatible with **Leica®, Thermo®, Tharmac®** and **Raymond Lamb®** staining machines.

AWS101-3-1M: Carries **Medite® 30 Place Racks**. Compatible with **Leica®, Thermo®, Tharmac®** and **Raymond Lamb®** staining machines

AWS101-3-2: Carries **Winlab® SR20, Sakura®, Hacker 20 Place Racks**. Compatible with **Leica®, Thermo®, Tharmac®,** and **Raymond Lamb®** staining machines.

AWS101-3-2D: Carries **InstrumeC®** and **Dako® 19 Place Racks**.

AWS101-3-2M: Carries **Medite® 20 Place Racks**. Compatible with **Leica®, Thermo®, Tharmac®** and **Raymond Lamb®** staining machines.

AWS101-3-9: Carries **InstrumeC®** and **Dako® 19 Place Racks**. Compatible with **Leica®,** and **Sakura®** dual staining

AWS101-3-17 Carries one only **InstrumeC®** and **Dako® 19 Place Racks**. Compatible with **Leica®,** and **Sakura®** dual staining




WINLAB PTY LTD

2 Pinnacle St, Brendale, Queensland, Australia, 4500

PO Box 5007, Brendale, Queensland, Australia, 4500

Ph: +61 7 3205 1209

Email: info@winlab.com.au www.winlab.com.au



Histology – Cytology – Adaptor Staining Hangers



AWS101-3-10: Carries **Winlab® SR20, Sakura®, Hacker®** 20 Place Racks. Compatible with **Leica®, and Sakura®** dual staining

AWS101-3-18: Carries one only **Winlab® SR20, Sakura®, Hacker®** 20 Place Racks. Compatible with **Leica®, and Sakura®** dual staining

AWS101-3-11: Dual Stainer compatible with **Winlab® SR30 Place Racks Dual**





AWS101-3-14: Carries **Winlab® SR20 Series Slide Racks**. Compatible with **Medite®** Staining machines

AWS101-3-16: Carries **InstrumeC® and Dako® 19 Place Racks**. Compatible with **Medite®** Staining machines

AWS101-3-5: Stainless Steel Handle Compatible with all 30 Place Slide Racks.

AWS101-3-6: Stainless Steel Handle Compatible with all 20 Place Slide Racks.



AWS101-3-12 SR 20 Plastic Handle; Compatible with all 20 Place Slide Racks.

AWS101-3-13 SR 30 Plastic Handle; Compatible with all 30 Place Slide Racks.