


WINLAB PTY LTD

2 Pinnacle St, Brendale, Queensland, Australia, 4500
PO Box 5007, Brendale, Queensland, Australia, 4500
Ph: +61 7 3205 1209
Email: info@winlab.com.au www.winlab.com.au



Winlab® Histology and Cytology Slide Racks

**OEM – Original and Replacement Slide Racks to match and integrate many standard slide-staining and cover-slipping machines.
Extremely rigid, and highly resistant to solvents up to 200° C.**



AWS101-2-1 SRLS30: A good all-rounder used by machine manufacturers; Vee Grooves for the hanger are vertical, horizontal and 45°; Compatible with **LEICA®**, **THARMAC®**, **RAYMOND LAMB®**, and **THERMO®** Staining Machines.
Max Temps 200°C



AWS101-2-2 SRM30: Exclusively made to suit the **MEDITE® Series** Staining Machines; Thinner body to fit narrow troughs; Includes the swivel feature found on those type of machine.
Max Temps 80°C



AWS101-2-3 SRLSM20: For machine integration, having a narrow bottom half and a larger top half, this rack is ideally suited to the narrow troughs in the **MEDITE®** Staining Machines and **SAKURA®** type Glass Cover-slipping machines.
Vee grooves for the hanger are vertical, horizontal and 45°.
Max Temps 200°C



AWS101-2-4 SRM20: Exclusively made to suit the **MEDITE®** Staining Machines, having a slim body to fit narrow troughs and includes the swivel feature found on **MEDITE®** Machines. Max Temps 200°C.



AWS101-2-5A SRLS20 (White): Ideally suited to the following cover slipping machines: **PHILIPS®**, **HACKER®**, **LEICA®**, **SAKURA®**, **THARMAC®**, **THERMO FISHER®**, **RAYMOND LAMB®**, **MEISEI®**, **VENTANA®**, **BAYER®** Staining and Cover-slipping Machines.

AWS101-2-5B (Black) SRLS20: Especially suited to Digital Pathology Slide Scanning Machines and most staining machines.
Vee Grooves for the hanger are vertical, horizontal and 45°.
Max Temps 200°C

AWS101-2-8 SRLM30: A good all-rounder used by machine manufacturers, having a swivel action hanger compatible with **MEDITE®**, **LEICA®**, **THARMAC®**, **RAYMOND LAMB®**, and **THERMO®** Staining Machines.
Max Temps 200°C

