

**WINLAB PTY LTD**

2 Pinnacle St, Brendale, Queensland, Australia, 4500  
PO Box 5007, Brendale, Queensland, Australia, 4500

Ph: +61 7 3205 1209

Email: [info@winlab.com.au](mailto:info@winlab.com.au) [www.winlab.com.au](http://www.winlab.com.au)




## Winlab® Large Slides Rack Insert

**NEW**

**Automated whole organ and multi-slide tissue staining on your regular machine with the Winlab Large Slides Rack Insert AWS101-3-3**


**Smooth edged design – Does not break or tear tissue from samples.  
Extremely durable, and highly resistant to solvents up to 190° C.**

For advice on your application contact [info@winlab.com.au](mailto:info@winlab.com.au)



### How to use:

- Simply place the **AWS101-3-3 Large Slides Rack Insert** into any rack at the desired spacing.
- Easily place the slides into the **Insert's** smooth track with no sharp edges to break the slides or tear your samples
- Large Slides Rack Insert (pair): Cat# AWS101-3-1



AWS101-2-5A – SR20 with Insert – all slides up to 2½" x 3"; 63mm x 76mm

AWS101-2-1 – SR30 with Insert – all slides up to 3" x 4"; 76mm x 101mm